# Vitique

## Esthetic Cementation System





CE 0482

090399/08.2011

### I. Vitique-System: Use of Base syringe

#### Use of base syringe: light cure of Vitique



Please refer to detailed description in the following instructions for use.




### II. Vitique-system: Assembly of catalyst syringe

#### Use of base + catalyst syringe: dual cure of Vitique



##### Removal of syringe caps

##### Assembly of base + catalyst syringe


##### Attaching of Smartmix-tip


#### Please note!

It is essential to check the fill level of the catalyst syringe. For this, you can use the ratcheting system on the piston rod. When the ratcheting system is not visible, the catalyst syringe is completely empty.

When the ratcheting system is not visible, the catalyst syringe is completely empty.

The ratcheting system is visible. The catalyst syringe contains material.

Please refer to detailed description in the following instructions for use.

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## Instructions for use

## English

Vitique is a highly esthetic radiopaque composite luting cement that can be used either light or dual cure. Vitique comes with try-in pastes that allow for the exact simulation of the final color of the restoration in the mouth.

In the light cure version, only the base is being used. In the dual cure version, a catalyst must be used in addition to the base. Vitique comes with a catalyst in two different viscosities, high viscosity and low viscosity. With different viscosities to choose from, there is the right viscosity available in Vitique for each given indication.

The mixing of the base and the catalyst in the Vitique system is accomplished by simply clipping the catalyst syringe to the base syringe (see also Vitique system, assembly of the catalyst syringe). The Smartmix-Tip allows for automatic mixing to an optimum quality and for direction application of the material.

#### Indications

Universal, permanent luting cement for highly esthetic cementation of ceramic or composite veneers, inlays, onlays, crowns or bridges; Vitique can also be used for the cementation of metal-based restorations.

#### Vitique-System

##### I. Use of base syringe - light cure only

- Remove cap (Fig. 01) and attach Veneer-Tip (Fig. 02 + 03). Please close syringe immediately after use to avoid premature setting of the material owing to ambient light (Fig. 04).

##### II. Assembly of catalyst syringe - dual cure

1. Lay out the base and catalyst syringe and a Smartmix-Tip, and remove caps from the syringes (Fig. 05 + 06)).

**Do not discard syringe caps since they will be needed later on for re-closure.**

2. Insert the base syringe at an angle from the top, with the outlet in the oval opening in the adapter plate at the outlet of the catalyst syringe (Fig. 07).
3. Press the grip plates of the basis and catalyst syringe together so as to ensure that the locking mechanism at the catalyst syringe snaps over the grip plate of the base syringe (Fig. 08).
4. Push the plunger of the catalyst syringe into the opening in the pressure plate of the plunger of the base syringe (Fig. 09).
5. Attach the Smartmix-Tip (Fig. 10). Ensure that the opening on the syringe is aligned with the opening on the tip.
6. Turn the tip 90° clockwise until it locks into position (Fig. 11).
7. Express a pea size amount of material and discard.

**The material is mixed in the tip upon extrusion and can be applied directly.**

**For disassembly of both syringes, follow steps in reverse order.**

#### Recommended use

1. Remove temporary restoration. Clean preparation with fluoride-free prophylaxis paste.

**Since the use of a rubber dam may change the tooth color by e.g. desiccation, it is recommended to determine the color of the luting cement with the try-in pastes before placing the rubber dam.**

2. Try in restoration.
3. Apply the transparent try-in paste onto the prepared tooth or the restoration. Seat restoration to place. If a color correction is desired, clean the restoration and the preparation with water and use the appropriate shade of the try-in paste to obtain the final desired color.

**The use of the catalyst paste for permanent cementation causes a slight color change in the base paste. Therefore also mix the try-in paste catalyst shade and the desired shade to simulate exactly the final cementation (see also: assembly of catalyst syringe).**

4. Remove the restoration. The restoration and the preparation are thoroughly rinsed with water spray.

**Since any residue of the try-in paste may affect the permanent cementation of the restoration, please ensure that the restoration and the preparation are thoroughly cleaned.**

**Ceramic surfaces that are silane treated by the laboratory become disturbed during try-in. Clean contaminated surfaces with alcohol or acetone and repeat application of silane.**

5. Etch internal surface of ceramic restorations with hydrofluoric acid. Please follow the instructions for use of the used hydrofluoric acid. Please note that some ceramics, like zirconium oxide, cannot be etched.
  6. Apply a drop of silane to the internal aspect of the restoration. We recommend the use of Silane (DMG). Please follow the instructions for use of the used silanating agent.
  7. Place rubber dam. If necessary, etch according to your preferred etching technique, e.g. Total-Etch.
  8. Apply the adhesive according to its instructions. We recommend the use of LuxaBond TE (DMG) for Vitique dual-curing and TECO (DMG) for Vitique light-curing.
  9. Apply Vitique to the internal aspect of the restoration and/or to the preparation (e.g. in case of an inlay cementation). When using the dual cure mode of Vitique, the working time is at least 2:00 minutes.
  10. Gently press restoration fully to place.
  11. Remove the bulk of the excess with a scaler or explorer. In case of veneer restorations, we recommend the use of a cotton roll or artist brush to remove excess cement with a gentle motion going from the restoration to gingival tissue.
- The excess material may also be exposed to a halogen lamp for approximately 2-3 seconds. This will slightly polymerize the excess material to allow for easier cleanup.**
12. Light cure the restoration for at least 40 seconds on each side of the tooth.
  13. Remove excesses, shape and finish with e.g. diamond burs, a #12 scalpel blade, finishing strips and flexible disks of increasing grits. High surface luster may be obtained with polishing paste using polishing cups/finishing disks.
  14. Evaluate the occlusion in the mouth and adjust as necessary.

#### A special note about the oxygen-inhibited layer

Upon polymerization of Vitique, areas that are exposed to oxygen produce an oxygen-inhibited layer. This layer does not set. The formation of this layer can be prevented as follows:

- Remove only the bulk of excess material first and leave a small excess. The oxygen-inhibited layer will then be removed together with the small excess when finishing the margins.

#### Please note

- Light curing units should have 450nm output and should be checked regularly. The light intensity of the units should be at least 400 mW/cm<sup>2</sup>. Place the light as close as possible to the material.

- The use of one-bottle bonding systems is contraindicated if Vitique is used without exposing it to light (e.g. when used in areas that do not allow any light penetration).
- Do not use Vitique together with materials that contain eugenol. We recommend the use of a eugenol-free cement such as TempoCem®NE (DMG) for temporary cementation.
- Please follow the instructions for use of the other materials that are used together with Vitique.

#### Interactions

Avoid materials containing eugenol, moisture and oily air; they can inhibit the polymerization of the material at the contact area.

#### Contraindications

The use of Vitique is contraindicated when a dry working area or the recommended application technique is not possible. Do not use directly on exposed pulp or if an allergy to any component is known.

#### Side effects

There have been no systemic side effects to date. In singular cases, hypersensitivity to components of the material cannot be excluded. Avoid contact of paste with skin. If accidental contact occurs, wash immediately with copious amounts of water and consult a physician.

#### Times

##### Vitique base

Removal of excess material expose to light for 2-3 seconds  
Final setting expose to light for 40 seconds

##### Vitique base + Vitique catalyst

Working time at least 2:00 min  
Removal of excess material expose to light for 2-3 seconds  
without exposing it to light, after 4:00 min  
Final setting expose to light for 40 seconds  
without exposing it to light, after 7:00 - 8:00 min

#### Composition

Vitique: Barium glass in a Bis-GMA based matrix of dental resins. Additives, catalysts, pigments. The catalyst paste contains also benzoyl peroxide. The filler volume is 67% by weight = 46 % by volume (0.02 - 4 µm).

Vitique Try-in: Silica gel and matched pigments in water-soluble polyglycol matrix

#### Classification

Vitique fulfills the requirements of ISO 4049:2000 (including color stability).

#### Storage

Do not store Vitique base and try-in pastes above 25 °C/77 °F. Store Vitique catalyst at 4 - 8 °C/39 - 46 °F. Do not use after expiry date.

#### Packaging

##### Intro Kit

8 x 1 syringe @ 6 g paste Vitique base (A1\*, A2.5\*, A4\*, B1\*, Bleach Light, Transparent, White, Pink); 1 syringe @ 2.1 g paste Vitique catalyst, High Viscosity  
1 syringe @ 1.9 g paste Vitique catalyst, Low Viscosity; 8 x 1 syringe @ 3.9 g paste Vitique-Try-in (A1\*, A2.5\*, A4\*, B1\*, Bleach Light, Transparent, White, Pink);  
2 syringes @ 1.3 g paste Vitique-Try-in, Catalyst Shade; 20 Smartmix-Tips, 100 Veneer-Tips REF 213013

##### Refills

##### Vitique base

1 syringe @ 6 g Paste, 10 Veneer-Tips  
A1\* REF 213038 A2.5\* REF 213039  
A4\* REF 213040 B1\* REF 213041  
Bleach Light REF 213042 Transparent REF 213043  
White REF 213044 Pink REF 213045

##### Vitique catalyst

1 syringe @ 2.1 g paste, High Viscosity; 10 Smartmix-Tips REF 213046  
1 syringe @ 1.9 g Paste, Low Viscosity; 10 Smartmix-Tips REF 213047

##### Vitique-Try-in

1 syringe @ 3.9 g paste, 10 Veneer-Tips  
A1\* REF 213048 A2.5\* REF 213049  
A4\* REF 213050 B1\* REF 213051  
Bleach Light REF 213052 Transparent REF 213053  
White REF 213054 Pink REF 213055  
1 syringe @ 1.3 g paste, 10 Smartmix-Tips  
Catalyst Shade REF 213056

##### Accessories

100 Veneer-Tips REF 213057  
50 Smartmix-Tips Brown 4:1 REF 213058

\* Shades according to Vita - Vita is a registered trademark of Vita Zahnfabrik H. Rauter GmbH

Federal Law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

Keep away from children! For dental use only!

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