# Silatec

091204/#2-06.2011

## Gebrauchsinformation

## Deutsch

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## Instructions for use

## English

Silatec, a pourable duplicating material based on A-silicones, is characterized by exceptional dimensional accuracy and precision. Silatec's balanced shore hardness, its high elasticity as well as its high tear strength guarantee an easy removal of casts. Exceptional precision in the reproduction of details and the material's excellent shape constancy produce best duplicating results. Molds made of Silatec can be stored and poured out even after months and can be used several times, if necessary.

### Indication

Fabrication of plaster or investment material duplicates

### Recommended Use

1. Put the master model into a standard duplicating mold.  
Do not wet the model with water. The model should have a normal room temperature of 18–25°C (64–77°F).
2. Shake both bottles of Silatec well before use and mix in a 1:1 ratio for at least 30 seconds until a homogeneous color is achieved. Depending on the size of the model, approximately 100 to 200g mixed material will be needed.  
When automatic mixing and dispensing units are used, shake material well and use up within 3 - 4 weeks. Alternatively, fill only the amount needed.
3. Pour the duplicating material in a continuous thin stream onto the dry model.  
The mixed material has a working time of 4 minutes.
4. Allow Silatec to set for 30 minutes at room temperature (18–25°C / 64–77°F).  
Remove the model by briefly blowing it out with compressed air and carefully removing it.
5. The mold can be poured several times with standard model cements or plasters. When doing so, wait 30 minutes before the next duplicating process. We recommend using special purpose hard plaster (ISO 6873, Type 3; e.g. Moldano, Bayer), ultra hard dental plaster for forms (ISO 6873, Type 4; e.g. Fujirock II, GC) or a heat resistant investment material (e.g. Stellavest, GC).

### Technical Data

Mixing ratio	1:1
Mixing time*	0:30 minutes
Working time* (incl. Mixing)	4:00 minutes
Setting time*	30:00 minutes
Strain in compression	≈ 16%
Compression set	≤ 0.25%
Detail reproduction	≤ 20 µm
Shore-A-hardness (1h)	Silatec ≈ 19

\* Times refer to normal room temperature (72°F / 23°C) and a normal relative humidity of air (50%). Higher temperatures shorten, lower temperatures increase those times.

### Composition

Vinylsiloxanes, Silicon dioxide, food pigments, additives, platinum catalyst

### Classification

ISO 14356:2003, Type 2: Duplicating Material

### Lagerung

Do not store above 25°C / 77°F. Do not use after expiry date.

### Packaging

#### Refill

1 bottle base @ 1 kg, 1 bottle catalyst @ 1 kg Silatec REF 909389

#### Economy Pack

1 bottle base @ 5 kg, 1 bottle catalyst @ 5 kg Silatec REF 909390

#### Super Saver Pack

1 bottle base @ 10 kg, 1 bottle catalyst @ 10 kg Silatec REF 909391

Federal Law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

Keep away from children! For dental use only!

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