# Supertec

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**Product description**

Supertec is a light-curing impression tray material used for the production of customized impression trays.

**Indication**

Production of customized impression trays for precision impression taking in the field of partial and full dentures, as well as implantology

**Basic safety instructions**

- For dental use only!
- Keep out of reach of children!

**Side effects**

Do not apply the material in cases of allergies to any of the components or in the event of contact allergies.

**Instructions for use**

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm². Place the light as close as possible to the material.

**Recommended use**

- Block out undercuts on the model with a suitable material, e.g. silicone impres-sion material.
- Mark the tray’s edge with a pencil and make a place-marker out of modeling wax.
- Remove protective cover from the prefabricated tray and adapt the tray to the model.
- Remove excess tray material along the pencil line with a modeling knife.
- Form a tray grip from the excess tray material and attach it to the model.
- If necessary, smooth material junctions with petroleum jelly.
- Cure the impression tray on the model for approx. 3 minutes with a suitable light unit, e.g. Kulzer® XS\* (UVA-VIS blended light/stroboscope).
- Note:** During curing it is possible to apply a light-curing varnish e.g. Dry Coat to the surface of the impression tray so that it is cured at the same time as the tray material. This creates a dry, shiny surface.
- If the impression tray is not treated with light-curing varnish, remove the dis-persion layer with alcohol or acetone. If light-curing varnish is used this is not necessary.
- Let the place holder wax boil out.
- Remove the impression tray from the model and cure the underside for approx. 3 minutes with a suitable light source e.g. Kulzer® XS\* (UVA-VIS mixed light/stroboscope).

**Note:** During curing it is possible to apply a light-curing varnish e.g. Dry Coat to the surface of the impression tray so that it is cured at the same time as the tray material. This creates a dry, shiny surface.

**Note:** The impression tray can be subsequently treated with corundum stones or fine carbide-fine cross-cut carbide burs. Wear suitable breathing protection when performing machining.

**Storage and shelf life**

- Store at room temperature (15 to 25 °C/59 to 77 °F)!
- Protect against exposure to direct sunlight!
- Do not use after the expiration date.

**Composition**

Glass filler and pyrogenic silicic acid in an UDMA-based matrix of dental resins.

**Technical data**

Transverse strength	> 60 MPa
Transverse modulus	> 4000 MPa

**Available product sizes**

50 trays upper jaw pink	REF 210410
50 trays lower jaw pink	REF 210411

\* Product names are trademarks or registered trademarks of the respective company.

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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