LuxaBond

Total Etch

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English

Product description

LuxaBond-Total Etch is a dual-cure bonding agent for enamel and dentin, suitable for use with dual-cure and self-cure composites, such as LuxaCore Z-Dual. It is necessary to etch enamel and dentin surfaces prior to application. Due to the dual-cure nature of the material, LuxaBond-Total Etch is particularly useful when bonding posts in areas where light can't penetrate.

Indications

Dual-cure bonding agent used in conjunction with the total-etch technique for core build-ups, and adhesion of resin cements for posts, crowns, bridges, inlays and onlays.

Contraindications

Do not apply the material in cases of allergies to any of the components, in particular benzoyl peroxide, or in the event of contact allergies.

Basic safety instructions

- For dental use only!
- Keep out of reach of children!
- Avoid contact with eyes, skin and oral mucosa. In the event of inadvertent contact with eyes, immediately rinse the eyes thoroughly with plenty of water and consult a physician if necessary. In the event of inadvertent skin contact, wash the affected area thoroughly with soap and water.

Side effects

- Isolated cases of contact allergies to components of the material cannot be excluded.

Interactions

Products containing eugenol, for example temporary cements or temporary root canal filling materials, may inhibit the polymerization.

- Products containing eugenol and other substances that inhibit polymerization must not be used together with the material.

Timing

Working time of the mixture of bond A + bond B (storage in a dark place)	5 min
Reaction time of the etching gel to enamel	20 to 60 s
Reaction time of the etching gel to dentin	max. 15 s
Pre-bond work-in time	15 s
Mixing time of bond A + bond B	5 s
Bond A + bond B work-in time	20 s
LuxaBond-Total Etch exposure time	min. 10 s

Instructions for use

- Light-curing units should have an output of 450 nm and should be checked regularly. The light intensity should be a minimum of 400 mW/cm². Place the light as close as possible to the material.
- Always use etching gel on enamel and dentin surfaces.
- In order to prevent premature polymerization due to ambient light, the material must be added to the mixing palette immediately prior to application.
- Ensure that the etched or bonded surface does not become contaminated prior to each subsequent working step.
- Do not remove the oxygen inhibition layer caused by the application because this is important for the bonding process with the material applied on top.
- Observe the manufacturer's instructions for all materials that are used with LuxaBond-Total Etch.

Recommended use

Use in the root canal, for cementing root posts, for cementing inlays

- Prepare the root canal or cavity in accordance with the general rules for adhesive procedures. We recommend the use of a rubber dam.
- Etch dentin and enamel with 37% phosphoric acid gel or similar. Start the etching process with the enamel and then expand this to include the dentin. The reaction time on the enamel should be 20 to 60 s. The reaction time on the dentin must not exceed 15 s.

Please note: When using LuxaBond-Total Etch it is essential to etch the affected enamel and dentin surfaces.

- Rinse with water for a minimum of 15 seconds.
- Remove general water excess. After the rinsing process, an even moist layer must remain on the dentin surface requiring treatment. If necessary, wet the dentin surface using a lightly moistened micro brush.
- Use a brush to work 1 to 2 drops of pre-bond into the etched tooth structure for 15 s. Remove any excess, e.g. with a paper tip. Blow gently using oil and water-free air. The surface should appear evenly moist.
- Mix 1 to 2 drops of bond A and bond B in a trough in the mixing palette at a ratio of 1:1 for approx. 5 s and work into the tooth structure for 20 s using a micro brush.

- Carefully remove bond excess from the preparation, e.g. with a paper tip, and then gently blow with oil and water-free air until the surface appears evenly moist.
- Apply the attaching cement in accordance with the manufacturer's instructions.

Please note: The working time of the attaching cement may reduce due to the contact with LuxaBond-Total Etch.

Please note: The use of LuxaCore Z-Dual is recommended for the cementation of root posts.

Please note: The use of Vitique is recommended for the cementation of inlays.

Core build-ups, cementation of crowns and bridges

Please note: Due to the fact that excess material may collect between the tooth and matrix (pooling effect), if possible use a matrix only after the bond has been applied.

- Prepare the cavity in accordance with the general rules for adhesive procedures. Bevel the enamel edges. We recommend the use of a rubber dam. In deep cavitated areas located close to the pulp, provide protection in the form of a suitable liner, e.g. by using a small quantity of calcium hydroxide and a thin layer of glass ionomer cement.
- Etch dentin and enamel with 37% phosphoric acid gel. Start the etching process with the enamel and then expand this to include the dentin. The reaction time on the enamel should be 20 to 60 s. The reaction time on the dentin must not exceed 15 s.

Please note: When using LuxaBond-Total Etch, it is essential to etch the affected enamel and dentin surfaces.

- Rinse with water for a minimum of 15 s.
- Remove general water excess. After the rinsing process, an even moist layer must remain on the dentine surface requiring treatment. If necessary, wet the dentin surface using a lightly moistened micro brush.
- Use a brush to work 1 to 2 drops of pre-bond into the etched tooth structure for 15 s. Remove any excess, e.g. with a paper tip. Blow gently using oil and water-free air. The surface should appear evenly moist.
- Mix 1 to 2 drops of bond A and bond B in a trough in the mixing palette at a ratio of 1:1 for approx. 5 s and work into the tooth structure for 20 seconds using a micro brush.
- Carefully remove any bond excess and then gently blow with oil and water-free air until the surface appears evenly moist.
- If required, light-cure LuxaBond-Total Etch for a minimum of 10 s using a suitable light curing unit.
- Apply core build-up material (e.g. LuxaCore Z-Dual) or resin cements (e.g. Vitique) in accordance with the manufacturer's instructions.

Please note: The working time of the resin cements or the core build-up material may reduce due to the contact with LuxaBond-Total Etch.

Please note: The use of Vitique is recommended for the cementation of onlays and crowns.

Composition

- Pre-bond: Ethanolic aryl sulfonate solution
- Bond A: Hydrophilic Bis-GMA based resin matrix, catalyst
- Bond B: Hydrophilic Bis-GMA based resin matrix, benzoyl peroxide

Storage

- Store LuxaBond-Total Etch in a refrigerator at 2–8 °C (35–46 °F).
- Close bottles immediately after use.
- Do not use the material past the stipulated shelf life or if the fluid assumes a gelatinous consistency.
- Ensure that the lids of bond bottles are not switched in order to avoid any cross contamination.

Available product sizes

Intro Kit	
1 bottle @ 5 ml pre-bond, 1 bottle @ 5 ml bond A, 1 bottle @ 5 ml bond B, 1 syringe @ 2 ml DMG etching gel, 25 Endobrushes red, 25 Endobrushes black, 3 Luer-Lock-Tips, 1 mixing palette	REF 212384

Accessories	
Endobrushes red	REF 212385
Endobrushes black	REF 212386
White Luer-Lock-Tips	REF 109942
DMG etching gel	REF 212873

Federal law (USA) restricts this device to sale by, or on the order of a dentist, or other practitioner licensed by law of the state in which he or she practices to use or order the use of this device.

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